

PIPETTOR ACCURACY CHECKS (BY VOLUME)

PERFORMED SEMI-ANNUALLY

Pipettor ID	Pipettor Volume	Date	Initials	Individual Deliveries (ml)										Average
				1	2	3	4	5	6	7	8	9	10	

Acceptable Ranges:

1 ml: 0.95 – 1.05 ml

Tolerance $\pm 5\%$

5 ml: 4.75 – 5.25 ml

Use Class A graduated cylinder

10 ml: 9.5 – 10.5 ml